

Can the stainless steel tubing cutter wheel type accommodate different size tubing - FAQ

Yes, the stainless-steel tubing cutter (cat# 46001-01) has a screw which is tightened to secure the tubing and can be adjusted for a range of different tubing diameters. It can be used for any tubing diameter in the range of 1/32 inch (0.79 mm) to 1/4 inch (6.4 mm).



The Cutting wheel has a thickness of 0.36 mm.

Click [HERE](#) fore ordering information for the tubing cutter and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com