



## What is the maximum pressure the Cogent Column Filters should be used with - FAQ

**The pre-column, Cogent Column Filters can be safely used with a maximum pressure of 5,000 psi.**

[Click here to view more information about pre-column filters for HPLC](#)



Printed from the Chrom Resource Center  
Copyright 2025, All Rights Apply  
**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900  
Fax: (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)