



Are AQ brand nylon syringe filters compatible with DMSO - FAQ

Yes, MicroSolv's AQ™ brand of syringe filters for dissolution testing and HPLC are made with a polypropylene housing which is DMSO compatible. Our nylon membranes, both 0.45um and 0.22um are fully compatible with solutions containing DMSO. Nylon is also compatible with esters, alcohols, weak acids and bases.

Click [HERE](#) for ordering information and pictures of AQ™ brand of syringe filters.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com