



For Normal Phase HPLC which syringe filter do you recommend  
- FAQ

**PTFE is the best membrane type for use with non polar, organic solvents required for Normal Phase Chromatography NP.** Click [HERE](#) for AQ™ brand PTFE syringe filter ordering information and pictures.

#### Attachment

**MicroSolv filters equivalency study pdf** [Download File](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)