


## Is the MicroSolv AQ brand nylon syringe filter compatible with formamide - FAQ


**Yes, our MicroSolv AQ™ nylon syringe filters are compatible with formamide.**

Formamide, dimethylformamide, and acetamide, are all compatible with our MicroSolv AQ™ brand nylon syringe filters.


*Acetamide:*



*Dimethylformamide:*



*Formamide:*



Click [HERE](#) for ordering information and pictures of AQ brand of syringe filters.

Printed from the Chrom Resource Center  
Copyright 2025, All Rights Apply  
**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900  
Fax: (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)