

## What is the pressure limit of the stainless steel filter sparging stones - FAQ

**All the sparging stone products are rated to a maximum pressure of 1,000 bar used under standard conditions.**

However, it is important to note that the tubing used and the connection of the tubing to the device may not perform at this pressure rating. The limiting pressure factor of a filter sparging stone is not for the filters / spargers, but for the connections.

Click [HERE](#) for filter sparging stone ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)