



## What is the pH range of 316 stainless steel laboratory tubing fittings and sparging stones - FAQ

The useful pH range of 316 stainless steel tubing, fittings and sparging stones is 1 - 14.

Click [HERE](#) for more information about stainless steel tubing

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)