



Are stainless steel tubing or fittings autoclavable - FAQ

Yes, many scientists have successfully used 316 stainless steel after autoclaving tubing and or fittings with gamma irradiation, heat or ethylene oxide sterilization methods.

Click [HERE](#) for stainless steel tubing ordering information

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com