

What is the definition of a fitting for HPLC - FAQ

The dictionary defines a fitting as “a small, often standardized part used to join other parts together properly.”

This is very similar to what plumbers refer as fittings. For the use of chromatographers, fittings consist of **nuts with internal and external threads** for both high pressure and low pressure connections. Also, unions, crosses and tees and adapters are considered fittings for both high and low pressure connections.

Click [HERE](#) to view our HPLC fittings

Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com