



## What are the critical dimensions of ferrules for low pressure fittings - FAQ

The "Bore Size" of the ferrules for low pressure fittings is critical for proper connection and for [Catalog# 45216-F-ETFE](#) or [45216-F-PTFE](#) it is ID  $1.70 \pm 0.05$  mm.

Designed for Use with 1/16th" OD low pressure tubing.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)