

## What is the useful temperature range of PEEK laboratory products - FAQ

The useful temperature range of PEEK related products for chromatography is up to 100°C (tubing) and up to 125°C (fittings).

PEEK is sometimes rated to 250°C but... It will not be as stable (lower burst pressure for tubing) as it is normally and some solvents may interact stronger than under RT with PEEK.

*PEEK has to be used with caution with higher temperatures.*

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)