



Which guard column do I use with the Partisil SAX Columns - FAQ

We suggest you use catalog # 66900-AGC (Anion Guard Cartridge (5 pack)), for use with Partisil™ HPLC columns.

You will also need a holder which is our catalog # 66900-WCS ,WCS guard cartridge holder, for use with Partisil HPLC columns.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com