



## What is the maximum temperature for Cogent Diamond Hydride HPLC columns - FAQ

**The maximum operating temperature for the Cogent Diamond Hydride™ HPLC column is 60° C.**

However, Higher Temperature Columns are available on request with High Temperature all Stainless Steel Hardware.



[Click HERE for Cogent Diamond Hydride Ordering Information](#)

Printed from the Chrom Resource Center  
Copyright 2025, All Rights Apply  
**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900  
Fax: (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)