



## Cogent amino HPLC columns designations for reversed or normal phase - Tech Information

Even though the Cogent NP Amino™ HPLC column is an excellent choice for non-hexane, normal phase chromatography, it is also very useful for highly acidic compounds using methanol or acetonitrile in reversed phase.

### IMPORTANT NOTE:

When the Cogent Amino™ column is labeled for normal phase, (denoted as **Cogent Amino NP column**) it is shipped in 90/10 heptane/ethyl acetate.

When the Cogent Amino™ column is labeled for reversed phase, (denoted as **Cogent Amino RP column**) it is shipped in methanol so one can easily go into reversed phase solvents.



Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)