



What is the maximum pressure of your guard columns - FAQ

The maximum pressure is the same as the analytical columns.

The standard 4um columns are rated up to 6000 psi but no more than 4500 psi is recommended during the run. The 2.0™ phases are suggested at no more than 9000 psi.

A PEEK column coupler will suffice unless you are using UHPLC, which is then not recommended.

Click [HERE](#) for ordering information and pictures



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com