



What are the metals content in Cogent HPS silica - FAQ

Metal ppm	Levels specification	Found
Na	<10	<1
Al	<10	<1
Fe	<10	<1
Ti	<10	<1
K	<10	None detected



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com