



## What are Cogent HPS HPLC Columns - FAQ

**Cogent HPS™ is a stationary phase (*both bonded and non bonded*) that is designed for applications that require high purity silica.**

Metal content in Cogent HPS™ is extremely low and this produces very good stability and “deactivates” the silica for use with basic analytes. A very smooth and uniform silica surface results for mechanical strength and stability.

Cogent HPS™ is available for HPLC with unmodified silica, C18, C8, Amino, Cyano or Phenyl stationary phases.

[Click HERE for More Information about Cogent HPS HPLC Columns](#)



Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)