



Do you offer a 0.5um titanium or PEEK frit for the PEEK column hardware - FAQ

In addition to the standard 5um stainless steel or PEEK frits, we have 0.5um titanium frits available but not in PEEK.

[PEEK Column Hardware Product Page](#)



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com