



## How many injections should I expect on the HSQ Kit HPLC column - FAQ

Thousands... if you do not drop it, overheat it, put too much pressure through it over 4,000 psi, subject it to *pressure shock* or inject anything on the column that is not part of the HSQ Kit™.

Click [HERE](#) for HSQ Kit ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)