

How do I validate the Excel spreadsheet that comes with the HSQ Kit - FAQ

There is no need for the end-user to perform any validation of HSQ Kit™ Excel® spreadsheet.

The document is fully validated by the manufacturer of the template. A certificate of validation can be provided upon request if needed (see below attached document). The scope of the validation encompasses the design of the spreadsheet, that its calculations are accurate and perform as intended.

This of course assumes that the end-user's input values themselves are accurate. The spreadsheet is also password-protected against altering of cell calculations and its validation demonstrates this. If the password is used to alter any cell, then the validation is considered void.

Click [HERE](#) for PQ & HSQ Kits™ ordering information and pictures

Attachments

Certificate of Validation PQ 7.07 pdf [Download File](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com