



The Gradient Visualization Solution components in the HSQ Kit - Tech Information

The HSQ Kit's Gradient Visualization Solution which is made under GMP conditions with calibrated balances follows a strict protocol has 1.0 ml/mL **Uracil** and **Caffeine** in 10% acetonitrile/water with 0.1% acetic acid as its components.

Click [HERE](#) for HSQ Kit™ ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com