



Composition of linearity solutions in the HSQ Kit - Tech Information

The HSQ Kit's™ linearity solutions are manufactured under GMP conditions with calibrated balances and following a strict protocol and comprise of **caffeine** dissolved in 90:10 **acetonitrile**:0.1% **acetic acid**.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com