



Differences between the linearity solutions in the HSQ qualification kit - Tech Information

They are all different and with precise concentrations of **caffeine** that allows you to easily and accurately determine the linearity of your HPLC UV detector.

L3 solution is the only linearity solution that has a **uracil** peak present. This peak is used as a pump flow stability marker and has no effect on caffeine or linearity tests.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com