



## How many tests will I get out of one HSQ Kit for HPLC qualification - FAQ

Depending on how careful you are with transferring the solutions, each kit will produce **up to 30 tests**. These tests can be on one instrument or on 30 different instruments.

Click [HERE](#) for PQ Kit™ and HSQ Kit™ ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)