

What is the difference between verifications validation and qualification - FAQ

In an LCGC magazine article by Drs. Ira Krull and Michael Swartz they reference the International Organization for Standardization (ISO) as follows:

“Validation is confirmation, through the provision of objective evidence that the requirements for a specific intended use or application have been fulfilled” whereas

Verification is confirmation, through the provision of objective evidence that specified requirements have been fulfilled.” They go on to say that if you use a compendia method, you verify and if you a non compendia or original method, you validate.

The term **Qualification** is now reserved for “instruments” whereas validation and verification are for methods and software.

For more details and further explanations see LCGC North America, Volume 23 Number 10 October 2005, Page 1100.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com