

Choosing a Proper Buffer in Capillary Electrophoresis of Proteins

The diverse types and characteristics of proteins makes selecting a proper buffer very challenging. Factors to consider when choosing a buffer for proteins are:

- A. The protein's pI and pK's
- B. The solubility of the protein
- C. The amount of glycosylation
- D. The effect pH has on the folding or unfolding of the protein

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com